English

SURFCATH™ REF 5590.106
6Fr – 0,8x2,0 mm – L. 200 mm

INDICATION AND DESCRIPTION

Surfcath™ is a surfactant injection catheter for the treatment of **preterm or full term** newborns at high risk of, or suffering from, respiratory distress syndrome; caused by lung surfactant deficiency (hyaline membrane disease).

Surfcath™ is a 20 cm long PEBA catheter, with a 6 Fr external diameter, and centimeter markings to guide the insertion. The distal end of the device (last 2cm) is a blunted black tip angled at 30° degrees. It is made of a 15 mm soft tip, and a 5mm extra-soft tip. The proximal end of Surfcath™ is Luer-Lock. The theoretical dead volume of Surfcath™ is 0,2 mL.

CONTRAINDICATIONS

Absolute Contraindications:
• Imminent need for intubation
• Maxillo-facial, tracheal or known pulmonary malformations
• Alternative cause for respiratory distress e.g. congenital pneumonia

Relative Contraindications:
• Severe RDS with high oxygen requirements, severe respiratory acidosis and/or widespread atelectasis on chest x-ray.
• Pneumothorax requiring drainage
• Prominent apnea, periodic or irregular breathing pattern

INSERTION

This procedure should only be performed by personnel thoroughly trained in this procedure or ET tube placement, and approved by the institution.

Prepare the necessary equipment:
- Laryngoscope
- Surfcath™ - Surfactant injection catheter
- Luer syringe pre-filled with surfactant. Always follow manufacturer's instructions in terms of storage and dosage of the surfactant.

Check the expiration date and the integrity of the packaging before starting the procedure. Observe sterile technique during the whole procedure.

1. Open the packaging and manually bend the Surfcath™ to the desired angle
2. Perform laryngoscopy using the appropriate size of blade
3. Insert the Surfcath™ orally or nasally with the concavity of the distal tip following the natural orientation of the respiratory airways.
4. Slide the Surfcath™ between the vocal cords using Magill forceps if necessary. Pass the epiglottis carefully thanks to the soft angled tip
5. The black intubation depth mark should be behind the vocal cords. Stop the progression of the catheter as soon as the black mark disappears.
6. Remove the laryngoscope while keeping the catheter in place. Close the mouth of the baby, and take note of the marking on the catheter which is at the lips or nostril of the baby. This will help you to know that the catheter stays in place.
7. Keep the catheter in place during the entire procedure.

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