


# Set up in three easy steps

## Mostcare Up


**Step one:** Connect Mostcare Up



Connect the Y cable between the bedside monitor and the ABP transducer

or

Connect the analogue cable between Mostcare Up and the bedside monitor




**Step two:** Enter patient data

1. Press **New monitoring**
2. Enter patient weight and height
3. Select ABP transducer and tick the **Y cable** box if applicable
4. Check other settings and press **Begin monitoring**



## Step three: Zero the transducer

- Place the transducer at the level of the phlebostatic axis
- Open the transducer stopcock to measure the atmospheric pressure
- Touch the curve on the screen and press **Zero ABP**
- Zero the primary monitor and re-establish trace



Haemodynamic variables		Physiological range**	Units
Pressures			
Dic	Dicrotic pressure	70 ÷ 105	mmHg
PP	Pulse pressure	30 ÷ 50	mmHg
MAP-Dic	Mean and dicrotic pressure difference	-10 ÷ +10	mmHg
Cardiac output			
SV	Stroke volume	60 ÷ 100	mL
SVI	Stroke volume index	35 ÷ 45	mL/m <sup>2</sup>
SVkg	Weighted stroke volume		mL/kg
CO	Cardiac output	4.0 ÷ 8.0	L/min
CI	Cardiac output index	2.6 ÷ 3.8	L (min· m <sup>2</sup> )
SVR	Systemic vascular resistance	800 ÷ 1400	dyne· sec/cm <sup>5</sup>
SVRI	Systemic vascular resistance index	1600 ÷ 2400	dyne· sec· m <sup>2</sup> /cm <sup>5</sup>
Oxygen delivery			
SpO <sub>2</sub>	Arterial oxygen saturation	96 ÷ 100	%
DO <sub>2</sub>	Oxygen delivery	900 ÷ 1000	mL/min
DO <sub>2</sub> I	Oxygen delivery index	500 ÷ 600	mL/min/m <sup>2</sup>
Efficiency and cardiac function			
dP/dt <sub>max</sub>	Arterial oxygen saturation	0.9 ÷ 1.3	mmHg/msec
CCE	Arterial oxygen saturation	-0.2 ÷ 0.3	units
CPO	Arterial oxygen saturation	0.80 ÷ 1.20	W
CPI	Arterial oxygen saturation	0.50 ÷ 0.70	W/m <sup>2</sup>
Vascular function			
Ea	Arterial elastance	1.10 ÷ 1.40	mmHg/mL
Dynamic variables			
SVV	Stroke volume variation	< 15*	%
PPV	Pulse pressure variation	< 15**	%
SPV	Systolic pressure variation	-	%

\* Approximate value reported in the literature in patients receiving controlled mechanical ventilation

\*\* Normal values in adult patients. Values will vary based on the clinical condition(s) that the patient has

Please note: For the most unstable patients, the number of beats considered for the average calculation should be increased: Touch the curve on the screen and increase the mean value.

 CRITICAL CARE

# Scan with your smartphone

to view online or visit [vygon.co.uk/mostcareevaluation/](http://vygon.co.uk/mostcareevaluation/)

Tel: 01793 748800 Copyright Vygon (UK) Ltd 2021 Content correct as of: 03/2021 Code: V00612 v3

