



Value Life

Neonatal
Complete Neonatal Nasogastric Feeding Range

Nutrisafe[®] 2

Supporting the delivery of
best practice

Nutrisafe[®] 2

Supporting the delivery of best practice

Nutrisafe[®] 2 is our complete and unique system that helps you meet the NPSA guidelines for safe, secure enteral delivery.



Safe

Nutrisafe 2 is incompatible with standard 6% luer taper devices. The lockable non-luer system reduces the risk of misconnection errors with no change to normal practice. The Nutrisafe 2 connector and the plastic cap tether have been designed to ensure it is safe for even your smallest patients.

Secure

The lockable system provides additional security. Accidental disconnection is virtually impossible. To disconnect the system you do not need to use force, you simply unscrew the device.

Confident

Colour-coded for easy recognition and size identification. The Nutrisafe 2 range features purple hubs and size related colour-coded caps, reducing the risk of accidentally attempting to connect devices intended for use with IV lines.

The unique Nutrisafe 2 secure lockable non-luer system gives you peace of mind that your patients are safe from misconnections and wrong route drug errors¹.

Comfortable


Designed with little babies in mind. The smooth shape was chosen specifically to eliminate sharp edges that could cause discomfort and markings on delicate skin.

30%

less dead space than the nearest competitor's device².

Think Small!

The compact design of Nutrisafe[®] 2 connections minimises dead space. This ensures safer, more accurate delivery of critical drugs such as caffeine, morphine, dexamethasone, chlorothiazide and spironolactone.



Think small! Nutrisafe 2 actual size compared to its nearest rival.



The Complete Solution

We offer a wide range of syringes, extension tubes, drawing-up straws, accessories and syringe pumps.

The Nutrisafe 2 syringe pump is an easy to use, high quality pump designed specifically for enteral feeding. The pump is capable of more than six hours of continuous work and is compatible with all types of Nutrisafe 2 syringes.



References

1. Patient Safety Alert 19. National Patient Safety Agency (NPSA). March 2007
2. Nutrisafe 2 Technical Report. Vygon (UK) Ltd, January 2013.

All references available on request.



www.carbonbalancedpaper.com

CBP00031881503134738

For further information, please contact: vygon@vygon.co.uk

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

Vygon (UK) Ltd • The Pierre Simonet Building • V Park • Gateway North
Latham Road • Swindon • Wiltshire • SN25 4DL

Tel: 01793 748800 • Fax: 01793 748899 • Twitter: [@vygonuk](https://twitter.com/vygonuk)

Web: www.vygon.co.uk

Content correct 03/2013 Code: 0062.NS2.NEO



www.vygon.co.uk