



Value Life

INFUSION SYSTEMS




**TIVA Sets
& A Sets**


PRODUCT GUIDE

TIVA Sets

Total intravenous anaesthesia from pump to patient




Features & benefits



Colour-coded Siamese tubing

For clearer identification of high and low flow drug administration.



Multi-lumen from hub to tip

Keeps infusions separate, reducing the risk of drugs mixing.



Armoured tubing


Eliminates the risk of kinking, ensuring continued infusions.

Ordering codes

| Vygon codes | NHSSC codes | Number of lumen | Gravity line | | High flow drug line | | Low flow drug line | | Box qty |
|-------------|-------------|-----------------|--------------|-------------------|---------------------|-------------------|--------------------|-------------------|---------|
| | | | Length (cm) | Priming vol. (mL) | Length (cm) | Priming vol. (mL) | Length (cm) | Priming vol. (mL) | |
| VAN7320BL | FSB1337 | Triple | 8.0 | 0.5 | 200 | 4.75 | 200 | 1.5 | 50 |
| VAN7330BL | FSB1338 | Triple | 8.0 | 0.5 | 300 | 7.0 | 300 | 2.5 | 50 |
| VAN7220BL | FSB1336 | Double | 8.0 | 0.5 | 200 | 4.75 | - | - | 50 |

A Sets

Total protection from pump to patient



Features & benefits



Anti-syphon valve

One-way valves that require a far greater driving pressure to open, minimising the risk of downloading.



Anti-reflux valve

Minimises the risk of back-tracking and enhances patient safety.



Halkey Roberts clamp

Secure, easy-to-use clamps, which enable one-handed operation during medication changes.

Ordering codes

| Vygon codes | NHSSC codes | Number of lumen | Detachable extension line | Gravity line | | High flow drug line | | Low flow drug line | | Box qty |
|-------------|-------------|-----------------|---------------------------|--------------|-------------------|---------------------|-------------------|--------------------|-------------------|---------|
| | | | | Length (cm) | Priming vol. (ml) | Length (cm) | Priming vol. (mL) | Length (cm) | Priming vol. (mL) | |
| VAN7300A | FSB1081 | Triple | No | 8.0 | 0.5 | 10 | 0.6 | 12 | 0.3 | 50 |
| VAN7300AL | FSB1082 | Triple | Yes | 8.0 | 0.5 | 185 | 5.0 | 187 | 4.7 | 50 |
| VAN7200A | FSB1079 | Double | No | 8.0 | 0.5 | 10 | 0.6 | - | - | 50 |
| VAN7200L | FSB1080 | Double | Yes | 8.0 | 0.5 | 210 | 5.1 | - | - | 50 |

For further information, please contact: vygon@vygon.co.uk

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

Vygon (UK) Ltd, The Pierre Simonet Building, V Park, Gateway North,
Latham Road, Swindon, Wiltshire SN25 4DL

Tel: 01793 748800 Fax: 01793 748899 Email: vygon@vygon.co.uk

vygon.co.uk [@vygonuk](https://twitter.com/vygonuk) [in](https://www.linkedin.com/company/vygonuk) vygonuk [v](https://www.youtube.com/channel/UCvYgONkUgQvYgONkUgQvYgONk) vygonuk [i](https://www.instagram.com/vygonuk) vygonuk

Copyright Vygon (UK) Ltd 2020 Content correct as of: 01/2020 Code: DXJB0005002.V00566 v1



Value Life