



Value Life

EMERGENCY



boussignac **CPAP+**[™]
Breathe easy, recover effectively

boussignacCPAP+™

Our trusted boussignacCPAP device + a new, improved mask for better fit and comfort

boussignacCPAP advantages

Safe

As an open system, boussignacCPAP is more comfortable for patients. The design ensures pressure in the lungs cannot be superior to the virtual valve so there is no risk of:

- Barotrauma/Volutrauma
- Hypoventilation

Easy to use

Straightforward regulation of FiO_2 with an integrated ring, if required. The flow of gas generates a pressure measured on the integral manometer.

A supplementary monitoring and administration port allows pressure measurement, EtCO_2 monitoring or additional oxygen or air delivery.

Effective

Patients who receive CPAP have improved oxygenation, lower respiratory rate, heart rate and systolic blood pressure.



Featuring...

A handy mask size guide on the outer packaging





Mask advantages

Created by Hans Rudolph Inc. a leading global designer and manufacturer of respiratory products, our boussignacCPAP+ offers the best emergency mask in the market, including:

Patient comfort

- Ergonomic and comfortable design for improved patient experience.

Innovation

- Innovative design providing lowest dead space vs alternative masks.

Oxygen retention

- Double seal, dual membrane mask, to prevent oxygen escaping.

Duration

- A true NIV mask that can be used for 24 hours.

Ordering information

boussignac**CPAP+** kit ordering codes



Standard kit

Our standard boussignac**CPAP+** kit contains our trusted boussignac**CPAP** plus our innovative new mask available in three sizes.



Mask size 4 (Small)
Vygon code:V08776615



Mask size 5 (Medium)
Vygon code:V08776606



Mask size 6 (Large)
Vygon code:V08776616

Pre-built kits

Our standard boussignacCPAP+ is also available in a range of pre-built kits for ease of use and ordering.



Kit option one

Our standard kit plus a FiO₂ regulator

Mask size	Code
S Small	V08776608
M Medium	V08776612
L Large	V08776614



Kit option two

Our standard kit plus a noise reducer

Mask size	Code
S Small	V08776605
M Medium	V08776613
L Large	V08776607



Kit option three

Our standard kit plus a FiO₂ regulator, noise reducer and nebuliser

Mask size	Code
S Small	V08776609
M Medium	V08776610
L Large	V08776611



Kit option four

Our standard kit plus a FiO₂ regulator, noise reducer, nebuliser and manometer

Mask size	Code
S Small	V08776617
M Medium	V08776618
L Large	V08776619

Ordering information

boussignacCPAP+ emergency kit ordering codes

boussignacCPAP+ emergency kits come in mini-cases for easy transportation and are suitable for all types of interventions. The kits include all the devices needed to administer CPAP to any patient.



Standard emergency kit

Our standard boussignacCPAP+ emergency kit contains our trusted boussignacCPAP plus four flowmeter options.

Configuration	Code
with Afnor flowmeter	5562.802
with DIN flowmeter	5562.842
with BSI flowmeter	5562.841
without flowmeter	5562.700

The boussignacCPAP+ emergency kits contain:



Manometer

Vygon code: 527.01



20ml syringe



Silicone harness

Vygon code: 5559.01



Mask size 4 (Small)

Vygon code: 5557.45



Mask size 5 (Medium)

Vygon code: 5557.55



Mask size 6 (Large)

Vygon code: 5557.65



Nebuliser

Vygon code: VBBNB31193



Elbow FiO₂ regulator

Vygon code: 5566.01



boussignacCPAP

Vygon code: 5570.13



O₂ flowmeter

Various devices available
(See page 7 for codes)

Accessories ordering codes

boussignacCPAP



boussignacCPAP
Vygon code: 5570.13



CPAP WEINMANN
Vygon code: 5570.132A

Single-use accessories



Manometer connector
Vygon code: 5558.053 or
Vygon code: 5558.203



E.T. tube connector
Vygon code: 555.01



vyconnector
Vygon code: 884.06



Noise reducer
Vygon code: 5558.91



Nebuliser
Vygon code:
VBBNB31193



Elbow FiO₂ regulator
Vygon code: 5566.01



Straight FiO₂ regulator
Vygon code: 5566.02



Star-shaped connecting tube
Vygon code: 5558.62

Reusable accessories



Afnor O₂ flowmeter
Vygon code: 5563.02



Afnor Air flowmeter
Vygon code: 5563.01



BSI O₂ flowmeter
Vygon code: 5563.41



BSI Air flowmeter
Vygon code: 5563.31



DIN O₂ flowmeter
Vygon code: 5563.42



DIN Air flowmeter
Vygon code: 5563.32



Nordic O₂ flowmeter
Vygon code: 5563.44



Nordic Air flowmeter
Vygon code: 5563.34



Manometer
Vygon code: 527.01



Standalone mini-case
Vygon code: 5562.799

Masks and harnesses



Mask size 3 (Paediatric)
Vygon code: 5557.35



Mask size 4 (Small)
Vygon code: 5557.45



Mask size 5 (Medium)
Vygon code: 5557.55



Mask size 6 (Large)
Vygon code: 5557.65



Fabric harness
Vygon code: 5559.03



Silicone harness
Vygon code: 5559.01

For further information, please contact: vygon@vygon.co.uk

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

Vygon (UK) Ltd, The Pierre Simonet Building, V Park, Gateway North,
Latham Road, Swindon, Wiltshire SN25 4DL

Tel: 01793 748800 Fax: 01793 748899 Email: vygon@vygon.co.uk

vygon.co.uk [@vygonuk](https://twitter.com/vygonuk) [in](https://www.linkedin.com/company/vygonuk) vygonuk [yt](https://www.youtube.com/channel/UCvYgONk8003V00447v2) vygonuk [ig](https://www.instagram.com/vygonuk) vygonuk

Copyright Vygon (UK) Ltd 2020 Content correct as of: 02/2020 Code: DXJB0008003 V00447 v2



Value Life