

ENTERAL

ENTERAL FEEDING PUMPS
AND ACCESSORIES

easymoov6

Intuitive nutrition



WHY CHOOSE VYGON?




**OVER
60
YEARS
EXPERIENCE**

**PATIENT CARE
LIES AT THE HEART
OF WHAT WE DO**



**CLINICAL
EXPERTISE**



WE ARE A
**CARBON
NEUTRAL
COMPANY**

MORE THAN
100 MILLION
PRODUCTS IN
110 COUNTRIES



CONTENTS

A new enteral feeding pump	3
An intuitive interface	4
A full range of giving sets	4
Blends project	5
Our commitment to the environment	6
References	6
Full product range and ordering codes	7

OUR SERVICE OFFERING



Customer Service & Technical Support

Talk through your enquiry with our dedicated Teams.



Education & Training

We offer our customers a variety of valuable and comprehensive training options to help you and your teams meet your training requirements.



A NEW ENTERAL FEEDING PUMP

FOR ALL PATIENTS ANYWHERE, ANYTIME

Designed to meet all patients needs in ambulatory and stationary medical care, the easymoov6 enteral feeding pump allows comfort and safety for patients and medical professionals in hospital and homecare.

Intuitive and available in different languages, this pump has a large screen showing all useful and relevant information.

Intuitive

All the information is available on one large screen: a short learning curve allowing the user to manage the pump quickly.

Quiet

A noiseless continuous rotary peristaltic rotor allows patients to feed during the night.

Powerful


A continuous rotary peristaltic rotor helps deliver high-viscosity food and preventing occlusion.

Safe

Nutrition modes, as well as visual and audible alarms, allow safe management of the pump by patients and/or caregivers.




AN INTUITIVE INTERFACE



A FULL RANGE OF GIVING SETS

The *easy*moovó pump has a full range of giving sets designed to cover all healthcare professionals' requirements and to improve safety and comfort thanks to the ENFit™ and nutri*safe*2 connections.

Some of the giving sets have an anti-free flow valve, allowing the hydration of patients.



BLENDS PROJECT

Blended ENteral Diet Study (BLENDS) is a project in conjunction with The University of Nottingham and a group of parents who use a blended diet. The Blends Project address the need for practical and user-friendly resources to support families, parents, and caregivers using a blended diet.

In November 2019, the British Dietetic Association (BDA) published a new, public facing, Policy Statement on the use of a blended diet with enteral feeding tubes which advised that a blended diet can be offered to tube fed individuals and their families as an option as part of patient centred and personalised care¹. However, there has yet to be much research to date that has aimed to establish the proportion of home enteral tube feeding (HETF) individuals in the UK using a blended diet.

BLENDS Project offers an interactive website and eBook that will guide parents of tube-fed children through their journey into the blended diet. The project has worked with medical experts and parents who currently use a blended diet, utilising a participatory co-creation approach to provide unique insight into the choices, challenges, and practicalities. This project fulfils a clear gap and will support clinical practice in the emerging speciality of blended diet.

THE NUMBER OF PEOPLE WHO ARE LONG-TERM, ENTERALLY TUBE-FED IN THEIR OWN HOMES HAS BEEN STEADILY INCREASING SINCE THE MID-1980S, AND BLENDING EVERYDAY FOODS WITH LIQUID TO BE GIVEN THROUGH AN ENTERAL FEEDING TUBE HAS BECOME INCREASINGLY POPULAR.



DISCOVER OUR JOURNEY
scan with your smart device



OUR COMMITMENT TO THE ENVIRONMENT

Vygon UK has been working towards creating a sustainable future for many years and we are a Carbon Neutral Company. We believe that sustainability isn't just about meeting current needs, but more importantly, it is about ensuring we are here for the long term and are paving the way for a bright tomorrow. We strive to deliver enhanced Corporate Social Responsibility (CSR), ensuring we manage the social, economic, and environmental effects of Vygon UK's operations responsibly in line with public expectations.



DISCOVER OUR JOURNEY
scan with your smart device



REFERENCES

- 1 British Dietetic Association. Policy Statement. The Use of Blended Diet with Enteral Feeding Tubes. 2019.

ORDERING CODES

Giving sets for easyMoov6 pump				
	Reference	Proximal connection	Type of giving set	Qty/box
Simple line giving sets	0VENP14S	ENPlus	ENFit™	30
	0VENP14S	ENPlus	ENFit™	30
	0VENP1GT	500 ml bag	ENFit™	30
	0VENP2GT	1200 ml bag	ENFit™	30
	0VENP3GT	1500 ml bag	ENFit™	30
Accessories	1MWPM6ACP	easyMoov6 AC adapter		1
	1MWPM6ACP	easyMoov6 Pole clamp		1
	0VEPM6BAG	Backpack for easyMoov6 pump		1
easyMoov6 pump	0VEPM6G02	easyMoov6 pump		1



09/2023 DXJ80002018 V01580 v1

FOR FURTHER INFORMATION, PLEASE CONTACT: info@vygon.co.uk

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature. This brochure has been printed on responsibly sourced and sustainable material. To help us reduce our carbon footprint all of our literature is available electronically either from your Product Specialist, on our website or by emailing info@vygon.co.uk

VYGON (UK) LTD, THE PIERRE SIMONET BUILDING, V PARK,
GATEWAY NORTH, LATHAM ROAD, SWINDON, WILTSHIRE SN25 4DL
RECEPTION: +44 (0)1793 748800 WWW.VYGON.CO.UK

 [@vygonuk](https://www.instagram.com/vygonuk)

