


Low Volume TIVA Set

Total intravenous anaesthesia from pump to patient



**NEW LOW
PRIMING
VOLUME**

Features & benefits




Anti-reflux valve

Minimises the risk of back-tracking and enhances patient safety.



Halkey Roberts clamp

Secure, easy-to-use clamps, which enable one-handed operation during medication changes.



Armoured tubing

Eliminates the risk of kinking, ensuring continued infusions.

Ordering codes

| Vygon codes | NHSSC codes | Number of lumen | Gravity line | | Low flow drug line | | Box qty |
|-------------|-------------|-----------------|--------------|------------------|--------------------|------------------|---------|
| | | | Length (cm) | Priming vol (mL) | Length (cm) | Priming vol (mL) | |
| VAN7220LV | FVR85306 | Single | 2.0 | 0.5 | 200 | 1.8 | 50 |